



CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL. ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538

AUTOCLAVES - 2 YEAR WARRANTY

FAST STERILISATION TIMES. DEEP CHAMBERS USING LESS DISTILLED WATER

18ltr or 22ltr CLASS B & S IN THE ONE UNIT

CHAMBER SIZE

18ltr - 250mm x 355mm

22ltr - 250mm x 450mm

European Class B with 3-time pre-vacuum preceding vacuum drying. It sterilizes and sanitizes miscellaneous packed or unpacked instruments whether solid or hollow. It efficiently rid air from all areas of any hollow instrument such as hand-piece making sure that water steam reaches every corner. It is installed with BOWIE & DICK test which measures the penetration of water steam. Printer attached to record sterilization process.



15ltr CLASS B



\$4500

CHAMBER SIZE - Diameter - 238mm x Depth - 370mm

European Class B with 3-time pre-vacuum preceding vacuum drying. It sterilizes and sanitizes miscellaneous packed or unpacked instruments whether solid or hollow. It efficiently rid air from all areas of any hollow instrument such as hand-piece making sure that water steam reaches every corner. It is installed with BOWIE & DICK test which measures the penetration of water steam. Printer attached to record sterilization process.

15ltr CLASS S




\$3999

CHAMBER SIZE - Diameter - 238mm x Depth - 370mm

A pre and post vacuum Class S autoclave that sterilizes solid instruments. It applies nozzle spray of water of water on chamber to generate steam fast and efficiently. Unit has fast sterilization times and a digital display. Printer is attached to record sterilization process.

12ltr CLASS B



\$3500

CHAMBER SIZE

Diameter - 230mm x Depth - 360mm

*** CLASS B INFO AS ABOVE ***



ADD GST TO ALL PRICING