

Capital Medical Supplies Pty Ltd

71 Heffernan Street, Mitchell. ACT. 2911

Ph. (02) 6241 5511

ABN 33 008 592 836

Fx. (02) 6241 7538

WELCH ALLYN EPISCOPE/DERMATOSCOPE



Fits all Welch Allyn 3.5v handles.

Features

- Precision, achromatic 10x optics system ensures crisp, clear lesion border definition and true pigment pattern rendition.
- Reticle contact plate facilitates measurement for more accurate diagnosis and objective tracking of lesion growth.
- Fits any 3.5v Welch Allyn power source for convenient, inexpensive additions to Welch Allyn diagnostic instrumentation.
- Large, knurled focusing ring allows you to get sharp focus easily.
- Rubber eyeguard provides superior comfort and folds down for spectacles.
- Owner's Manual includes introduction to skin surface/epiluminescence microscopy, tables of diagnosis criteria, and 50 detailed lesion photographs.