

# Capital Medical Supplies Pty Ltd

71 Heffernan Street, Mitchell. ACT. 2911

Ph. (02) 6241 5511

ABN 33 008 592 836

Fx. (02) 6241 7538

## "AURORA" 3 Metre Vision Chart

- Solid-State Illumination Panel with a life expectancy over 5000 hours
- Remote Control of Illumination for each Acuity Level
- Battery or Mains Operation
- Letter & Illiterate Charts included
- Enables Testing of Both Adults and Children
- Portable
- Very Cost Effective
- Patented **Electroluminescent** Technology



The "Aurora" system can be connected to the mains power using the transformer provided or can function with the four AA size batteries which are already installed in the battery tray.

