



**CAPITAL MEDICAL SUPPLIES PTY LTD**

**71 HEFFERNAN STREET, MITCHELL. ACT. 2911**

**Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538**

***WALL MOUNTED HEIGHT  
MEASURE***

***0 TO 200CM***



**Handy Height Scale**

**Handy Height Scale (PE027)**

**Ideal for taking consistent height measurements.**

**Fix the top portion of the unit to a wall 2 metres above the floor, then lower the bottom portion to the patient's head and read the measurement.**

**Range, 0-2000mm.**