



CAPITAL MEDICAL SUPPLIES PTY LTD
71 HEFFERNAN STREET
MITCHELL. ACT. 2911

PH. 02 62415511
FAX. 02 62417538
E-MAIL. info@capitalmedicalsupplies.com.au
ABN. 33 008 592 836

HEINE ILLUMINATED SIGMOID/PROCTO/ANOSCOPE SYSTEM WITH DISPOSABLE SCOPES



- Annular Fiber Optic illumination. Bright, reflex-free illumination.**
 - Magnifying lens 1,5x. High-resolution image.**
 - Hinged viewing window with insufflation port. Improved vision by insufflation.**
 - Stainless Steel.**
- Instrument head for UniSpec disposable tubes TAKES ALL 3 TYPES OF DISPOSABLE SCOPES**