

HI·bebe Fetal Doppler BT-200



Specification

Display

LCD

Heart Rate Range

Low Battery Display

Accuracy

Signal Quality Display

Ultrasound Frequency

50~240 bpm

Ultrasound Intensity

±2%

Speaker Output Power

2MHz

<10mW/cm²

1.2W

Audio output for earphone

Built-in Loudspeaker

Auto Shut-off

5 minutes

Auto Shut-off on signal

1 minute

Power

1.5V Battery X2 (Type:AA)

Sensitivity

16 Weeks Onward

PC Interface

Sound Card

Battery Life

360 minutes