



**CAPITAL MEDICAL SUPPLIES PTY LTD**  
71 HEFFERNAN STREET  
MITCHELL. ACT. 2911

**PH. 02 62415511**  
**FAX. 02 62417538**  
**E-MAIL. [info@capitalmedicalsupplies.com.au](mailto:info@capitalmedicalsupplies.com.au)**  
**ABN. 33 008 592 836**



## **Histofreezer—Low Cost Cryotherapy**

### **Portable Cryosurgical System**

**A patented cyrogenic liquid mixture of compressed gas comprised of dimethyl ether and propane. This gas mixture does not harm the ozone layer. The cryogen is released into a foam applicator cooling it to an effective freezing temperature of -57 degrees C.**

**A truly portable cryogenic system. May be left in travel bag for nursing home or other situations in which it is not convenient to take liquid nitrogen.**