



**CAPITAL MEDICAL SUPPLIES PTY LTD**

**71 HEFFERNAN STREET, MITCHELL. ACT. 2911**

**Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538**

## *Headlight Ri-Focus LED*



**With Lithium Batteries**



### **ri-focus® LED**

The headlights ri-focus® and ri-focus® reverse impress with their extraordinary suitability for everyday use. The bright and non-glaring illumination makes these devices both patient and user friendly to the highest degree. The lighting position of the ri-focus® is above the forehead and between the eyes for the ri-focus® reverse, in order to be suitable for all kinds of examination requirements.

Halogen lamp 8 V for white and even light.

The lens at the front of the light head can be focused to allow exact adjustment of the size of the illuminated examination area.

Integrated heat sink for patient-friendly use.

Step switch on cable for regulating light intensity in two steps.

Electrical cable runs along the outside of the headband to permit unhindered use.

Practical and fast exchange of light heads using special fixture located on the side of the forehead.

Simple adjustment and fixation of the headband with help of the rotating aluminium knob.

Soft, exchangeable foam protection around the forehead for comfort in wearing the headband.