

CMS**Capital Medical Supplies Pty Ltd**
71 Heffernan Street, Mitchell. ACT. 2911**CMS****Ph. (02) 6241 5511****ABN 33 008 592 836****Fx. (02) 6241 7538**



Littmann Model 3000 Electronic Stethoscope

Product Information

- Four times (400%) the amplification (15-20 dB) of the best conventional stethoscopes makes indistinct sounds easier to hear
- Filter circuitry provides enhanced clarity and three frequency response modes for optimal heart and lung auscultation:
Bell (20-200 Hz),
Diaphragm (100-500 Hz),
Extended Range (100-1,000 Hz)
- Easy-to-use ergonomically designed chestpiece. Virtually as light as the 3M™Littmann™ Master Cardiology stethoscope
- Excellent acoustic seal and comfortable fit with patented 3M™Littmann™ Snap Tight soft-sealing eartips
- Comfortably angled, anatomically correct headset
- Patient-friendly, nonchill rim
- One-year warranty
- Volume control provides stepless adjustment to preferred levels
- 40 hours of continuous use from one AAA alkaline battery
- Alloy Chestpiece sized for adult or pediatric patients