

Capital Medical Supplies Pty Ltd

71 Heffernan Street, Mitchell. ACT. 2911

Ph. (02) 6241 5511

ABN 33 008 592 836

Fx. (02) 6241 7538

RESCUAID ADULT RESUSITATOR SET



Features

Adult resuscitator 1500ml with 60cm H2O pop-off valve

2500ml Oxygen reservoir

Two durable adult masks #5145, #5144

Two guedel airways

See through carry case with instruction manual

Code	Description
------	-------------

REUSABLE ADULT MANUAL RESUSITATOR COMPLETE IN CASE.

ALSO AVAILABLE WITHOUT THE CARRY CASE - LESS EXPENSIVE

ALSO AVAILABLE AS DISPOSABLE RESUSITATOR
IN SIZES INFANT - CHILD - ADULT ALL \$49 EACH