

Capital Medical Supplies Pty Ltd

71 Heffernan Street, Mitchell. ACT. 2911

Ph. (02) 6241 5511

ABN 33 008 592 836

Fx. (02) 6241 7538

Proport Model 1600 Electronic Baby Scale



- 1 year guarantee
- Displays weight in both kilograms/grams and pounds/ounces
- Capacity 20kg or 44lb
- Graduation: 0 - 10kg : 10g; 0 – 20lb : 1/2oz; 11 – 20kg: 20g; 21 - 44lb: 1oz
- Large power saving LCD readout, 19.2mm (0.77”) digits
- Lock-in weight feature – a stable weight can be obtained even if the baby is moving
- **Tare weigh function**
- Removable baby tray for easy storage or weighing of a small child
- Auto zero & shut-off; low battery indicator
- Requires one 9 volt alkaline main battery (not included)

