



# CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538

*SECA HEAVY DUTY FLOOR SCALE 6" LARGE DIAL  
HOSPITAL USE  
MADE IN GERMANY*



Large Clear Display  
Ideal for remote, Power free areas  
Very Robust construction  
Non slip foot-tread  
Hardwearing and longlasting  
Very Economical

Capacity : 150kg  
Readability : 500g  
Power: N/A  
Colour: Black / White  
Weight: kg  
Dimensions: 300 x 105 x 460 mm