




CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL. ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538



seca

**MODEL 700 SLIDING WEIGHT
SCALE COMPLETE WITH
INTEGRAL HEIGHT MEASURE.**

THE MEASUREMENT CENTRE.

Professional Mechanical Scales MODEL 700

- 220kg Capacity !!
- Graduations 50gm
- Complete with height rod
- Height rod graduations 60-200cm
- Dimensions 438W x 1438H x 538Dmm
- Weight 19kg
- 12 Months Warranty
- Manufactured Germany
- Fantastic mechanical scale