

Capital Medical Supplies Pty Ltd – 71 Heffernan St, Mitchell, ACT. 2911

Ph. (02) 6241 5511 Fax. (02) 6241 7538

Welch Allyn® Spot Vital Signs® -- NIBP/Nellcor Pulse Oximetry/Temperature Model 42NTB-E1

Vital signs in seconds, payback in no time!

Taking vital signs just got easier with the Welch Allyn Spot Vital Signs, the first ever hospital-grade, automated, multi-parameter device created specifically to spot-check vital signs. The Spot Vital Signs offers measurements that are quick, easy, affordable, and accurate; lightweight portability; and a variety of mounting options to meet your needs.



Features

- Standard non-invasive blood pressure (NIBP) as well as pulse rate and MAP.
- Optional SureTemp® thermometry.
- Optional Nellcor's motion-tolerant SpO2 technology with finger clip sensor.
- Uses the oscillometric method to determine systolic and diastolic pressure, MAP and pulse rate.
- Adult durable one-piece blood pressure cuff.
- Large LCD display with backlight.
- Extremely user-friendly and easy to use.
- Lightweight and portable:

• Available on a Mobile Stand

Technical Specifications

Dimensions: 9.70"(24.64cm) x 5.72" (14.53cm) x 4.73" (12.01cm)

Weight: 2.04 kg without stand

Blood Pressure Accuracy: Meets or exceeds SP1--1992 AAMI Standards.

Warranty

2-year warranty on all parts and labor; 1-year warranty on temperature probe and SpO2 sensor.