



CAPITAL MEDICAL SUPPLIES PTY LTD
71 HEFFERNAN STREET
MITCHELL. ACT. 2911

PH. 02 62415511
FAX. 02 62417538
E-MAIL. info@capitalmedicalsupplies.com.au
ABN. 33 008 592 836



Batteries Included

- 1 second measurement
- Multiple Use (Wide Range Measurement)
- Accurate and reliable
- Gentle and Easy to Use
- Auto-Display Memory
- Safe and Hygienic
- Fever Alarm

Microlife Digital Infrared Ear Thermometer

Measuring Range: 0 °C to 100.0 °C (32.0 °F to 212.0 °F)

Accuracy: Laboratory: ±0.2 °C, 32.0 ~ 42.2 °C

Display: Liquid Crystal Display with indicating unit 0.1 °C

Acoustic: a. The unit is turned ON and ready for the measurement: 1 short «bi» sound.

b. Complete the measurement: 1 long beep sound.

c. System error or malfunction: 3 short «bi» sounds.

Memory: a. Auto-Display the last measured temperature when switched on. b. 12 readings recall in the Memory Mode.

Operating

temperature: 5°C to 40 °C (41.0 °F to 104 °F)

Storage/transport

temperature: -25 °C to +55 °C (-13 °F to 131 °F)

Automatic Switch-off: Foolproof Switch-off: Auto-off if START button is pressed and held for 15 seconds.

Battery: CR2032 BATTERY (X1) - at least 1000 measurements

Dimensions: 153mm (L) x 31 mm (W) x 40 mm (H)

Weight: 53g (w/ battery), 50g (w/o battery)

Standards: Complies with ASTM E-1965 requirements