



CAPITAL MEDICAL SUPPLIES PTY LTD
71 HEFFERNAN STREET
MITCHELL. ACT. 2911

PH. 02 62415511
FAX. 02 62417538
E-MAIL. info@capitalmedicalsupplies.com.au
ABN. 33 008 592 836



Features

- **3.5v PanOptic™ Ophthalmoscope with Axial PointSource Optics offers a 5X larger field of view and 25% more magnification vs. traditional ophthalmoscopes. Patented Glare Extinguishment System enables glare-free views of the retina.**
 - **The 3.5v MacroView(TM) Otoscope with built-in throat illuminator provides 30% more magnification and twice the field of view of traditional otoscopes, along with the ability to adjust the focus.**
 - **Lithium Ion handle provides twice the on-time with just half the weight of conventional handles for comfortable use. The low battery indicator insures you don't get caught with a "dead" handle.**
 - **Handle accepts all Welch Allyn 3.5V heads**

