



**CAPITAL MEDICAL SUPPLIES PTY LTD**  
71 HEFFERNAN STREET  
MITCHELL. ACT. 2911

**PH. 02 62415511**  
**FAX. 02 62417538**  
**E-MAIL. [info@capitalmedicalsupplies.com.au](mailto:info@capitalmedicalsupplies.com.au)**  
**ABN. 33 008 592 836**

## **Welch Allyn Pocket Oscopes & Sets**



### **Features**

- 2.5V Otoscope also has built in throat illuminator - two in one instrument! Features halogen illumination and fiber optics for a cool light with no reflections or obstructions, plus a rectangular viewing lens for instrumentation under magnification.
- 2.5V PocketScope handle has a smooth rheostat to control light intensity and features a handy pocket clip.

- **2.5V PocketScope Ophthalmoscope uses halogen illumination for true tissue color and consistent, long lasting illumination. It has 6 aperture selections for general and specialist use, and 48 lenses for better resolution.**
- **2.5V Otoscope also has built in throat illuminator - two in one instrument! Features halogen illumination and fiber optics for a cool light with no reflections or obstructions, plus a rectangular viewing lens for instrumentation under magnification.**
- **2.5V PocketScope handle has a smooth rheostat to control light intensity and features a handy pocket clip.**

