



CAPITAL MEDICAL SUPPLIES PTY LTD
71 HEFFERNAN STREET
MITCHELL. ACT. 2911

PH. 02 62415511
FAX. 02 62417538
E-MAIL. info@capitalmedicalsupplies.com.au
ABN. 33 008 592 836



Welch Allyn Wall Diagnostic Set.
Your choice of Heads.
Fits all 3.5v attachments.

- **Handle-based rheostats for quick light-intensity adjustments.**
- **Comfortable, contoured handles that are easy to hold and maneuver.**
- **Durable, easy to clean, high-impact ABS housing.**
- **Improved cradle design features OptiSense™, a reliable optical sensor that automatically turns handles on when they are removed from the cradles.**
- **Halogen HPX™ lamp provides 30% more light output.**
- **You may add a SureTemp® thermometer; KleenSpec® Plus Dispenser; third handle module for additional instrument; and a wall aneroid.**