

RESUSCITATION DEMAND VALVES & MASKS

OXYGEN DEMAND VALVE RESUSCITATOR

SPECIFICATIONS

Manual Flowrate: 40 – 60 l/min	Supply Pressure: 400 kPa nominal	Inlet Connection: 9/16" – 18 D.I.S.S.
Demand Flowrate: 0 – 160 l/min	Pressure Relief Valve: 54 – 60 cmH ₂ O	Patient Connection: 22 /15 mm male / female taper
Demand Flow Triggering System: -5 cmH ₂ O @ 100 l/min, maximum	Exhalation Resistance: -5 cmH ₂ O @ 50 l/min, maximum	Weight: 230 grams

NOTE: D.I.S.S. (Diameter Indexed Safety System as per CGA V-5)

- ◆ Safe and effective means of delivering oxygen
- ◆ Can be used as a demand valve resuscitator or as a standard demand valve
- ◆ Integral pressure relief valve
- ◆ Can be attached to any standard 400 kPa oxygen outlet, including: hospital outlet valves, regulators, or Twin-0-Vacs

Ordering Information:

Oxygen Demand Valve Resuscitator

Part No. 552095

Applications:

The Demand Valve Resuscitator provides a safe and effective means of delivering an oxygen concentration to ventilate non-breathing or respiratory distressed patients. The unit can also be used as a standard Demand Valve to deliver oxygen to spontaneously breathing patients upon demand. The unit is simple, compact and always ready for instant use and attaches to any standard 400 kPa oxygen outlet.

Spare Parts:

Spare Diaphragm Valve	Part No. 552082
Spare Patient Valve Assembly	Part No. 552083
Oxygen Hose Assembly, 2 mtr D.I.S.S – S.I.S connections	Part No. 552086



552095



552086



552083



552082

RESUSCITATION MASKS

These adult and child masks feature an integral cuff and shell and are designed to help you in emergency situations when spontaneous breathing has ceased.

They have a soft cushion and clear shell giving you both an effective seal and easy visibility of the patient's face and lips.

Ordering Information:

Resuscitation Mask, Adult, Size 5
Resuscitation Mask, Child, Size 3

Part No. 515743
Part No. 515744



515743



515744