

MULTEPOINT TERMINAL WALL OUTLETS

Part No.	Gas	Colour	Connection Mode	Outlet Connection
553040	Oxygen	White	Screw type / Single action	S.I.S Sleeve
553041	Medical Air	Black	Screw type / Single action	S.I.S Sleeve
553042	Nitrous Oxide	Blue	Screw type / Single action	S.I.S Sleeve
553042	Direct Suction	Yellow	Screw type / Single action	S.I.S Sleeve
553044	Venturi Suction	Yellow	Screw type / Single action	S.I.S Sleeve
553045	Direct Scavenge	Beige	Screw type / Single action	S.I.S Sleeve
553046	Venturi Scavenge	Beige	Screw type / Single action	S.I.S Sleeve
553047	Surgical Tools	Aqua	Screw type / Single action	S.I.S Sleeve
553048	Entonox	Blue/White	Screw type / Single action	S.I.S Sleeve
518900	Oxygen	White	Screw & Quick connect type / Dual action	S.I.S Sleeve
518901	Medical Air	Black	Screw & Quick connect type / Dual action	S.I.S Sleeve
518902	Nitrous Oxide	Blue	Screw & Quick connect type / Dual action	S.I.S Sleeve
518903	Direct Suction	Yellow	Screw & Quick connect type / Dual action	S.I.S Sleeve

NOTE: S.I.S. (Sleeve Indexed Safety System as per AS2896)

- ◆ High quality performance and reliability
- ◆ Available in a range of gases and suction services
- ◆ Unique gas specific diameter indexing system prevents incorrect assembly of units on installation or when dismantled for maintenance
- ◆ High flow capability
- ◆ Quick connect or screw type models available
- ◆ Minimum parts for simplicity and lower maintenance costs
- ◆ Highly visible colour coded name plates make it easy to identify the gas or suction service
- ◆ Suits various installation requirements – wall, pendant, column or duct mounting
- ◆ Gas supply tube can be positioned either up or down to suit installation (MK 3 models only)

Applications:

The Multepoint 3 ranges of terminal wall outlets are designed for hospital gas systems. This series of model continues the tradition for high quality, performance and reliability of the Multepoint name.

Accessories:

Oxygen Handwheel Kit	Part No. 518406
Air Handwheel Kit	Part No. 518405
Nitrous Oxide Handwheel Kit	Part No. 518407



ABOVE:
Multepoint Single Action



LEFT:
Multepoint Dual Action