

R-VIVA EMERGENCY RESUSCITATOR

- ◆ Simple, lightweight and resilient
- ◆ Plenty of extra storage pockets to store and protect your systems components
- ◆ Basic kit design allows for easy addition of other accessories or components
- ◆ Oxygen Demand Valve Resuscitator provides a resuscitation flowrate of 40-60 l/min and a demand flowrate of 0-160 l/min
- ◆ Suction unit provides a static vacuum of -300mmHg, with a 30 l/min maximum developed suction. With a 140ml trap jar and a nipple to suit 8mm suction tubing
- ◆ Pin indexed yoke (Fig.2) gas inlet connection
- ◆ Therapy flow attachment offering a fixed flowrate of 8 l/min or 15 l/min with a nipple to suit 5mm therapy tubing
- ◆ R-Viva bag made of Fabrec 1000 Denier material
- ◆ Weighs approximately 4kg without a cylinder fitted (approximately 7kg with cylinder fitted)
- ◆ Dimensions are 64 x 28 x 25cm

Ordering Information:

R-Viva Resuscitator with
Demand valve & Suction unit

Part No. 552160

Applications:

The R-Viva Emergency Resuscitator is designed for the professional rescuer. The main component of the kit is the Multi-Flow piston regulator. This regulator has two oxygen outlets to run both the Comweld Demand Valve Resuscitator and Suction unit as well as oxygen therapy selector with two flow settings. This means that the complete kit can provide benefits of resuscitation, suction and oxygen therapy from the one compact unit. Operating off a 400 litre (C-size) oxygen cylinder, these three different treatment modes can be used simultaneously or separately for almost any emergency situation.

Spare Parts:

Resuscitation Mask, Adult, Size 5	Part No. 515743
Resuscitation Mask, Child, Size 3	Part No. 515744
Oxygen Demand Valve Resuscitator	Part No. 552095
Oxygen Supply Hose for Demand Valve Resuscitator or Suction unit	Part No. 552086
Comweld Suction Unit	Part No. 554052
Spare Bodok Seals for cylinder (pk of 10)	Part No. 551022



552160
(Cylinder not included)

RESUSCITATION RESUSCITATION KITS

R-VIVA OXYGEN THERAPY KIT

- ◆ Simple, lightweight and resilient
- ◆ Comfortable carry handle and shoulder strap
- ◆ Stretcher attachment straps to secure the unit during transportation
- ◆ Internal pocket space
- ◆ External velcro strips for the attachment of an additional pouch, if necessary
- ◆ Non slip cylinder holder
- ◆ Clear viewing window allows a direct view of the cylinder contents gauge and flow setting when in use
- ◆ Therapy flow control offering ten flow selections at flowrates of 1,2,3,4,5,6,7,8,10 and 15 l/min and a nipple to suit 5mm therapy tubing
- ◆ R-Viva bag made of Fabrec 1000 Denier material
- ◆ Weighs approximately 1.5kg without a cylinder fitted (approximately 4.5kg with cylinder fitted)
- ◆ Dimensions are 55 x 28 x 15cm

Ordering Information:

R-Viva Oxygen Therapy Kit

Part No. 552162

Applications:

The R-Viva Oxygen Therapy kit is ideal for use in private residences, workplaces and elderly citizens homes for the supply of oxygen. The main component of the kit is the Multi-Flow piston regulator, which has an oxygen therapy selector with ten flow settings. Operating off a 400 litre (C-size) oxygen cylinder this kit is great for almost any oxygen therapy requirement. The inclusion of a mouth to mask system for emergency resuscitation situation adds to this units versatility.

Spare Parts:

Resuscitation Mask, Adult, Size 5	Part No. 515743
Resuscitation Mask, Child, Size 3	Part No. 515744
Spare Bodok Seals for cylinder	Part No. 551022



552162
(Cylinder not included)