

SUCTION SUCTION ATTACHMENTS

COMWELD SUCTION ATTACHMENTS - Central Vacuum

Part No	Colour	Gauge Range	Connections	
			Inlet	Outlet
554006	Yellow	0 to -100 kPa	S.I.S Vacuum Handwheel	Nozzle to suit 8mm I.D. Tubing
554007	Yellow	N/A	S.I.S Vacuum Handwheel	Nozzle to suit 8mm I.D. Tubing

NOTE: S.I.S (Sleeve Indexed System Handwheel as per AS2896)

- ◆ Designed to suit a Central Vacuum system
- ◆ Horizontal positioning of connector and gauge to the suction nozzle, for wall mounting
- ◆ On / Off control tap
- ◆ Suction nozzle and filter are fitted as standard and are also available as spare parts
- ◆ Easy to read dual scale (kPa and mmHg) vacuum gauge (applicable models)
- ◆ Colour coded pressure gauges and identification labels for quick recognition of the gas / function of the product
- ◆ Gas specific inlet fittings to prevent fitting to the wrong gas service
- ◆ Suction outlet connection is a nominal 8mm I.D. tubing nipple
- ◆ Gauge range is 0 to -100 kPa
- ◆ Weighs approximately 500 grams

Applications:

These Comweld Suction Attachments are a basic connector with no controls other than an On/Off tap. They simply serve as the adaptor between the suction performance of the pipeline outlet / central suction system and the receiver jar assemblies (or other suction collection apparatus being used).

Spare Parts:

Spare Sintered Filter	Part No. WA185
Spare Outlet Nozzle/Filter Holder	Part No. 518634
Spare Vacuum Gauge (0 to -100 kPa)	Part No. 522501
Spare O-ring (for outlet nozzle)	Part No. 7964024



WA185



518634



SUCTION SUCTION ATTACHMENTS

COMWELD SUCTION ATTACHMENTS - Venturi Suction

Part No	Colour	Gauge Range	Connections	
			Inlet	Outlet
554013	Yellow	0 to -100 kPa	S.I.S Vacuum Handwheel	Nozzle to suit 8mm I.D. Tubing
554010	Yellow	N/A	S.I.S Vacuum Handwheel	Nozzle to suit 8mm I.D. Tubing

NOTE: S.I.S (Sleeve Indexed System Handwheel as per AS2896)

- ◆ Designed to suit a Venturi Suction system
- ◆ Non-return valve
- ◆ Horizontal positioning of connector and gauge to the suction nozzle, for wall mounting
- ◆ Suction nozzle and filter are fitted as standard and are also available as spare parts
- ◆ Easy to read dual scale (kPa and mmHg) vacuum gauge (applicable models)
- ◆ Colour coded pressure gauges and identification labels for quick recognition of the gas / function of the product
- ◆ Gas specific inlet fittings to prevent fitting to the wrong gas service
- ◆ Suction outlet connection is a nominal 8mm I.D. tubing nipple
- ◆ Gauge range is 0 to -100 kPa
- ◆ Weighs approximately 400 grams

Applications:

These Comweld Suction Attachments are a basic connector with no controls. They simply serve as the adaptor between the suction performance of the venturi pipeline outlet and the receiver jar assemblies (or other suction collection apparatus being used).

Spare Parts:

Spare Sintered Filter	Part No. WA185
Spare Outlet Nozzle/Filter Holder	Part No. 518634
Spare Vacuum Gauge (0 to -100 kPa)	Part No. 522501
Spare O-ring (for outlet nozzle)	Part No. 7964024



554013



554010

COMWELD LOW SUCTION ATTACHMENTS - Venturi Suction

Part No	Colour	Gauge Range	Connections	
			Inlet	Outlet
518535	Yellow	0 to -25 kPa	S.I.S Vacuum Handwheel	Nozzle to suit 8mm I.D. Tubing

NOTE: S.I.S (Sleeve Indexed System Handwheel as per AS2896)

- ◆ Designed to suit a Venturi Suction system
- ◆ Safety valve designed to prevent negative pressures greater than -25 kPa being supplied
- ◆ Easy to read dual scale (kPa and mmHg) vacuum gauge which gives clear indication of low pressures
- ◆ Non Return Valve which prevents any positive pressure being applied should a fault occur in the suction source
- ◆ A positive pressure relief valve within the patient circuit for additional safety
- ◆ Horizontal positioning of connector and gauge to the suction nozzle, for wall mounting
- ◆ Suction nozzle and filter are fitted as standard and are also available as spare parts
- ◆ Colour coded pressure gauges and identification labels for quick recognition of the gas / function of the product
- ◆ Gas specific inlet fittings to prevent fitting to the wrong gas service
- ◆ Suction outlet connection is a nominal 8mm I.D. tubing nipple
- ◆ Gauge range is 0 to -25 kPa
- ◆ Weighs approximately 500 grams

Applications:

This Comweld Suction Attachments is a basic connector with no controls. It serves as a low suction adaptor between the venturi pipeline outlet and the receiver jar assemblies (or other suction collection apparatus being used). This unit is designed to keep the suction performance of the pipeline outlet to a lower level than normal, for more specific low suction applications.

Spare Parts:

Spare Sintered Filter	Part No. WA185
Spare Outlet Nozzle/Filter Holder	Part No. 518634
Spare O-ring (for outlet nozzle)	Part No. 7964024

