

ACCESSORIES & SPARES

3 WAY ADAPTORS



Extends the use of a single gas outlet to enable two gas hoses or two pieces of equipment to be used simultaneously. All connections coded and self sealing.

Ordering Information:

Oxygen 3 Way Adaptor
Air 3 Way Adaptor
Nitrous Oxide 3 Way Adaptor

Part No. 551140
Part No. 551141
Part No. 551142

HOSE ASSEMBLIES



Ready to use gas supply hoses consisting of 4 metres of tubing fitted with colour coded, non interchangeable handwheels.

Ordering Information:

Oxygen Hose Assembly
Air Hose Assembly
Nitrous Oxide Hose Assembly

Part No. MA164
Part No. MA165
Part No. MA166

HANDWHEEL KITS

Each kit contains a colour coded non-interchangeable handwheel plus three different sized nipples (5mm, 1/4", 8mm), O-rings, O-clips and plastic shrink sleeve.

Ordering Information:

Oxygen Handwheel Kit
Air Handwheel Kit
Nitrous Oxide Handwheel Kit

Part No. 518406
Part No. 518405
Part No. 518407



GAUGE



Ordering Information:

Contents Gauge, Plain, 30000 kPa, Series-O Regulator
Contents Gauge, Nitrous Oxide, 20000 kPa, Series-O Regulator
Contents Gauge, Air, 30000 kPa, Series-O Regulator
Contents Gauge, Oxygen, 30000 kPa, Series-O Regulator
Flow Gauge, Vacuum, -100kPa, Twin-O-Vac & Suction Attachments
Flow Gauge, Oxygen, 0-14 l/min, TM17 Flowmeter

Part No. 551042
Part No. 551043
Part No. 551044
Part No. 551045
Part No. 552501
Part No. 551041

522501

MISCELLANEOUS

ACCESSORIES & SPARES

CYLINDER KEYWHEEL OR SPANNER



511956



551162

A keywheel or spanner design product to turn medical gas cylinders on and off with a minimum of effort.

Ordering Information:

Cylinder Keywheel
Cylinder Spanner

Part No. 511956
Part No. 551162

FLOWMETER NUT / NIPPLES (Packet of 10)



515830



515849

Suits Ezi-Flow and TM style flowmeters

Ordering Information:

Oxygen Nut / Nipples (pack of 10) **Part No. 515830**
Air Nut / Nipples (pack of 10) **Part No. 515849**

BODOK SEALS (Packet of 10)



551022

Suits Series-0 regulators or any equipment with pin indexed yoke inlet fittings for medical gas cylinders.

Ordering Information:

Bodok Seals (Packet of 10)

Part No. 551022

BACTERIA FILTERS (Packet of 100)



554050

Spare filters for all Twin-0-Vac units.

Ordering Information:

Bacteria Filters (Packet of 100) **Part No. 554050**