



**CAPITAL MEDICAL SUPPLIES PTY LTD**  
71 HEFFERNAN STREET  
MITCHELL. ACT. 2911

PH. 02 62415511  
FAX. 02 62417538  
E-MAIL. [info@capitalmedicalsupplies.com.au](mailto:info@capitalmedicalsupplies.com.au)  
ABN. 33 008 592 836

**UNION  
2 SECTION  
ELECTRIC  
HI-LOW  
COUCH**



- \* Grey base / Navy upholstery / 710mm wide
- \* Two sections: 680mm [+80°] / 1220mm
- \* Electric height adjustment [45 - 95cm] with foot control
- \* Two leg / swing arm direct lifting system with bearings on moving parts.
- \* Single action retractable castors.
- \* Lifting capacity 200kg
- \* Matches RACGP standard.

