




CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL, ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538



UNION ELECTRIC HI-LOW COUCH

- ***Grey base / TBA upholstery / 600mm wide****
 - * ***Two sections: 680mm / 1220mm (+80°)***
 - * ***Electric height adjustment (50 - 100cm) with foot control***
 - * ***Two leg / swing arm direct lifting system with bearings on moving parts.***
 - * ***Single action retractable castors.***
 - * ***Lifting capacity 175kg***
- Delivery charges apply***