




# CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL. ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538



- \* Grey base / TBA upholstery / 600mm wide\*
- \* Three sections: 360mm with breathe hole & plug (+40 / -90deg) / 560mm / 940mm (+80deg).
- \* Electric height adjustment (50 - 100cm) with foot control
- \* Two leg / swing arm direct lifting system with bearings on moving parts.
- \* Single action retractable castors.
- \* Lifting capacity 175kg



### *The Union 3 Section Chair/Couch*

*This couch is very versatile as it can be used for the patient to be sitting up or lying flat.*

*Delivery is normally 4-6 weeks from time of order.*

*Delivery Charges apply*