




CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL, ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538



- * Grey base / TBA upholstery / 710MM wide*
- * Three sections: 360mm with breathe hole & plug (+40 / -90deg) / 560mm / 940mm (+80deg).
- * Electric height adjustment (50 - 100cm) with foot control
- * Two leg / swing arm direct lifting system with bearings on moving parts.
- * Single action retractable castors.
- * Lifting capacity 200 kg



The Union 3 Section Chair/Couch

This couch is very versatile as it can be used for the patient to be sitting up or lying flat.

Delivery is normally 4-6 weeks from time of order.

Delivery Charges apply

**MEETS NEW RACGP
GUIDELINES**