



CAPITAL MEDICAL SUPPLIES PTY LTD
71 HEFFERNAN STREET
MITCHELL. ACT. 2911

PH. 02 62415511
FAX. 02 62417538
E-MAIL. info@capitalmedicalsupplies.com.au
ABN. 33 008 592 836



Liquid Withdrawal Device

The Liquid Withdrawal Device enables the LD25, LD35 or the LD50 to conveniently pressure transfer liquid nitrogen at rates up to 8 liters per minute. The low [10 psig (0.7 bar/69kPa)] assures safe dispensing through a phase separator or an accessory transfer hose. A smaller phase separator is available.



Roller Base & Tipping Stand

The roller base is a sturdy, aluminum, five-wheel dolly that enables movement of refrigerators and dewars easily and safely. Choose from three sizes to fit either the 18.8 in (478 mm) dia., 15.6 in. (396 mm) dia., or the 26.7 in (678 mm) dia. units. The tipping provides mobility and assist in transferring liquid nitrogen for the LD25 dewar.