



CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL. ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538

Jamar Hydraulic Hand Dynamometer

Gives accurate grip strength readings without "feeling" the handle move.

The adjustable handle has 5 positions to accommodate any size hand. Registers up to 91kg and has dual scales that show kg and pounds.



Advanced Hand Dynamometer

Advanced Hand Dynamometer (PE022)

This dynamometer is an upgraded version of the PE007.

The frame is stronger and sturdier with the hand grip and general operation being smoother and more user friendly.

Comes complete with sturdy case and instructions. Measures a grip strength of up to 100kg.

