





**CAPITAL MEDICAL SUPPLIES PTY LTD**  
**71 HEFFERNAN STREET**  
**MITCHELL. ACT. 2911**

**PH. 02 62415511**  
**FAX. 02 62417538**  
**E-MAIL. info@capitalmedicalsupplies.com.au**  
**ABN. 33 008 592 836**

USE AT 3 METRES (Eye Chart)



**Snellins Eye chart Standard.**  
**Available as 3 Meter or 6 Meter Test**  
**British Optical—Wall Hung Type**  
**Non Illuminated.**



**Snellins Eye chart Illiterate.**  
**Available as 3 Meter or 6 Meter Test**  
**British Optical—Wall Hung Type**  
**Non Illuminated.**

**Charts are available for very young children or illiterate adults that do not require letter recognition. One version uses simple pictures.**