



# CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL, ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538



The Hyfreicator<sup>®</sup> has been the #1 product of choice by physicians in office based electrocautery for over 60 years. CONMED is pleased to continue this tradition with the new Hyfreicator<sup>®</sup>2000 which features:

- Memory to store up to 3 power settings
- Remote control handpiece to change power settings, maintain the sterile field and decrease procedural time
- Dual Microprocessors: one to control the function and output of the Hyfreicator<sup>®</sup>, the second to terminate activation if a hazard or safety issue is detected
- Cost effective alternative to various laser procedures
- Available with height adjustable stand
  
- The Hyfreicator 2000 comes complete with
- Sharp & blunt tips
- Non-sterile & sterile pencil sheaths
- Up/down pencil; instructional video
- Wall mounting kit

When it comes to your patients care . . . accept nothing short of Excellence