






**CAPITAL MEDICAL SUPPLIES PTY LTD**  
71 HEFFERNAN STREET  
MITCHELL. ACT. 2911

PH. 02 62415511  
FAX. 02 62417538  
E-MAIL. [info@capitalmedicalsupplies.com.au](mailto:info@capitalmedicalsupplies.com.au)  
ABN. 33 008 592 836

## Ishihara Color Test Plates



Ishihara's test is composed of a series of colored cards on which numbers or lines of equal shade can be read by a person with normal color vision but not by someone with defective color vision. It makes use of the peculiarity that in red-green blindness, blue and yellow appear remarkably bright compared to red and green. The color plates are encased in specially designed album-type book for ease of handling. It includes four special plates for tests to determine the kind and degree of defect in color blindness. This color blindness test is accepted by leading authorities.

Available in Book form of 14 Plate, 24 Plate, 36 Plate

