

Oscilla[®] SM910 Screening Audiometer



*** Pulse tone ***
*** Easy to operate ***

SPECIAL FEATURES

Pulse tone.
Very readable displays.
Noise free push buttons.
11 frequencies: 125 Hz - 8000 Hz.
Calibration without entering the cabinet.

Oscilla® SM910 Screening audiometer is a modern manual audiometer, which you may learn to operate within a few seconds.

Design as well as advanced technology has been enforced, and by means of processor controlled technique, the **Oscilla® SM910 Screening audiometer** is still made for the future. For the tone generator is used technology, which assures a 100% stable tone without any delay or distortion.

Aluminium cabinet provides a neat product, which at the same time is very sturdy and protects in the best way against any environmental influences from other electronic equipment; as aluminium protects completely against environmental hazards.

Technical specifications:

Frequencies: 11 frequencies from 125 Hz to 8000 Hz.
Hearing Level Range: -10 dB to 90 dB in 5 dB increments.

Maximum intensities:

Frequencies Hz	125	250	500	750	1000	1500	2000	3000	4000	6000	8000
Air conduction dB	50	70	90	90	90	90	90	90	90	80	70

Program: Manual test.
Calibration: Calibration data stored in programmable memory. (EN60645-1).
Power supply: Netadapter 12V AC (220-240V).
Dimensions: 225mm x 180mm x 55mm.
Weight: 650g. (unit only).
Standard accessories: TDH-39 earphones, SILENTA noise reducing headset, patient response, netadapter.
Optional extras: Carrying bag.
