



# CAPITAL MEDICAL SUPPLIES PTY LTD

71 HEFFERNAN STREET, MITCHELL. ACT. 2911

Ph. (02) 6241 5511 ABN. 33 008 592 836 Fax (02) 6241 7538

## SCHILLER-Switzerland *CARDIOVIT AT-2*

Fully automatic, microprocessor controlled 6-channel A4 electrocardiograph, mains and battery operated.

### Standard unit includes:

- 12 simultaneous standard leads
- Pacemaker detection
- High resolution thermal printer
- Direct function keyboard with rubber keys
- Printout on A4 recording paper
- Battery with life for up to 400 recordings between charges

### Standard Accessories

Including Software "C" Interpretive Software

- ECG measurements
- Average complexes and measurement markers
- Detailed measurement results table
- Interpretation statements for children and adults
- 10 lead patient cable with banana plugs
- 500 Bio-adhesive electrodes SCHILLER-Tabs Ag
- Snap Clip adaptors
- Recording paper
- Mains cable



PRINTOUT NOT TO SCALE—NORMAL IS A4 SIZE

